

## Chitosan membranes for spatially controlled cell adhesion and specific cell recruitment

C. A. Custódio<sup>1,2,3</sup>, A. M. Frias<sup>1,2</sup>, A. del Campo<sup>3</sup>, R. L. Reis<sup>1,2</sup>, J. F. Mano<sup>1,2</sup>

<sup>1</sup>3B's Research Group - Biomaterials, Biodegradables and Biomimetics, University of Minho, AvePark, Zona Industrial da Gandra, S. Cláudio do Barco, 4806-909 Caldas das Taipas - Guimarães, Portugal.

<sup>2</sup>ICVS/3B's, PT Government Associated Laboratory, Braga/Guimarães, Portugal.

<sup>3</sup>Max-Planck-Institut für Polymerforschung, Ackermannweg 10, 55128 Mainz, Germany.

We propose a concept of biomaterials that are able to fix specific cell types onto their surface when in contact with a mix population of cells.

Adipose tissue has shown to be an interesting source of stem cells with therapeutic potential. However only a small amount of the heterogeneous mixture of the cells extracted from lipoaspirates are stem cells, and within stem cells there are different populations with different capabilities to differentiate through a lineage.

We studied the ability of immobilized antibodies on chitosan surfaces to capture specific types of cells with a spatial micrometer resolution.

Antibodies were covalently immobilized onto chitosan membranes using bis[sulfosuccinimidyl] suberate (BS3). X-ray photoelectron spectroscopy (XPS) was used to chemically characterize the surface and quartz crystal microbalance (QCM) to calculate the amount of adsorbed and/or immobilized antibody.

Data shown greater immobilization when BS3 was used compared to simple adsorption. Specific antibodies covalently immobilized in a surface, kept their bioactivity and controlled the type of cell that attached on the chitosan surface. Microcontact printing permitted to covalently immobilize antibodies in patterns allowing a spatial control in cell attachment. Cell sorting experiments performed using a mixture of adipose stem cells and osteoblast like cells shown that chitosan surfaces were able to capture a specific phenotype depending on the immobilized antibody.

Selected Symposium: Specific Biomaterials as Basis of TERM, Natural Derived Polysaccharides